

## Specifications

EndoStratus™ Irrigation Pump

REF#: EGA-500

Dimensions: 200(W) x 173(H) x 385(D)mm

Weight: 11.4 lbs

Power Requirements: 100-240 V

Max Flow Rate:

Auxiliary water port: 0-280 ml/min

Biopsy channel: 0-600 ml/min

EndoStratus™ CO<sub>2</sub> Insufflator

REF#: EGA-501


Dimensions: 200(W) x 173(H) x 385(D)mm

Weight: 12.5 lbs


Power Requirements: 100-240 V

## Easy Set Up and Installation

① Choose a CO<sub>2</sub> Source




Wall Regulator



CO<sub>2</sub> Tank

② Connect Hose and Tubing to Insufflator



③ Ready for Procedure

## Compatible Accessories

### Tubing



REF #100130  
EndoGator® Tubing

REF # 100145CO2  
Endo SmartCap® Tubing for Olympus®  
140/160/180 Series Endoscopes

REF # 100165CO2  
Endo SmartCap® Tubing for  
Fujinor® Endoscopes and Pentax®  
Endoscopes

REF # 100150CO2  
Endo SmartCap® Tubing for Olympus®  
100/130 Series Endoscopes and  
Pentax® Endoscopes



REF #100603  
EndoGator® Kit Tubing

REF #100611  
EndoGator® Kit Tubing with CO<sub>2</sub>

REF #100609  
EndoGator® Hybrid Tubing with CO<sub>2</sub>

Compatible with Olympus® Endoscopes

### Connectors



REF #100241  
Single Use Auxiliary Water Jet Connector for Olympus® Endoscopes

REF #100242  
Single Use Auxiliary Water Jet Connector for Pentax® Endoscopes

REF #100141  
Endo SmartCap® Adapter for Fujinor® 530 Series Colonoscopes

### Biopsy Irrigation



REF #100303  
DEFENDO® Y-Opsy® Irrigator  
for Olympus® or Fujinor® Endoscopes

MEDIVATORS®

ENDO  
**stratus™**

IRRIGATION PUMP and CO<sub>2</sub> INSUFFLATOR



MEDIVATORS®, EndoStratus™, EndoGator®, Endo SmartCap®, DEFENDO®, and Y-Opsy® are trademarks or registered trademarks of Medivators Inc.  
Fujinor® is a registered trademark of Fuji Photo Film Co., Ltd. and/or its affiliated entities.  
Olympus® is a registered trademark of Olympus Corporation, Olympus America, Inc., and/or their affiliated entities.  
PENTAX® is a registered trademark of PENTAX of America, Inc.

**MEDIVATORS®**  
A Cantel Medical Company

MEDIVATORS Headquarters  
14605 28th Avenue North  
Minneapolis, MN  
Tel: +1.800.444.4729  
Fax: +1.763.553.3387

MEDIVATORS BV  
Sourethweg 11  
6422 PC Heerlen  
The Netherlands  
Tel: +31.45.5.471.471  
Fax: +31.45.5.429.695

MEDIVATORS Asia/Pacific Pte Ltd  
1 A International Business Park  
#05-01 Singapore 609933  
Tel: +65.6227.9698  
Fax: +65.6225.6848

MEDIVATORS Beijing Representative Office  
Room 708, 7th Floor Kaiheng Center, Block B  
No. 2 Chaoyangmennei Street  
Dongcheng District, Beijing China 100010  
Tel: +8610.6567.8446  
Fax: +8610.6567.8445

[www.medivators.com](http://www.medivators.com)

**MEDIVATORS®**  
A Cantel Medical Company

## The CO<sub>2</sub> Solution for GI Endoscopy

The MEDIVATORS® EndoStratus™ CO<sub>2</sub> Insufflator is a versatile system that connects to either a wall or tank CO<sub>2</sub> source. Combined with the full line of the disposable MEDIVATORS® Endo SmartCap® CO<sub>2</sub> Tubing, the EndoStratus™ CO<sub>2</sub> Insufflator offered by MEDIVATORS® allows clinicians to offer the optimal procedure for the patient.

MEDIVATORS® offers a complete range of CO<sub>2</sub> solutions to suit individual facility needs for Olympus®, Pentax® and Fujinon® GI endoscopes. With direct connection to sterile water and any CO<sub>2</sub> tank or wall source - MEDIVATORS® has the complete solution for your facility and patients.

- CO<sub>2</sub> is easier on the patient both during and after the procedure due to rapid absorption by surrounding tissue.<sup>1</sup>
- Minimizes patient discomfort associated with bloating and cramping from room air.<sup>1</sup>



Daily Use Endo SmartCap® CO<sub>2</sub> Source Tubing

## Warm Water Technology


Deliver irrigation for GI endoscopy procedures with a warm water option by utilizing the EndoStratus™ Irrigation Pump. Combine the pump with any GI endoscope and irrigate through the auxiliary water jet or biopsy channel to improve procedural visibility.

- Greater visibility for improved diagnosis and polyp detection
- Auxiliary water jet or biopsy channel irrigation capability
- Warm water technology
- Compatible with MEDIVATORS® EndoGator® Tubing and Connectors
- Easy to use water bottle cradle
- Automatic prime button provides instant irrigation upon foot pedal depression



Auxiliary Water Jet Connection

## CO<sub>2</sub> Insufflator Features



## Irrigation Pump Features

